

BEL PHOTONICS ► DV-5000



BEL Photonics DV-5000 camera is a High-resolution color camera studied for those applications where high performance is required. It is provided with high speed USB2.0 interface and large area CMOS image sensor, that gives a large field of view. These features together with included user-friendly Image analysis software BEL Eurisko allow capture of high quality images with excellent colors and resolution.

DV-5000 è un modello di telecamera ad alta risoluzione appositamente studiata per le applicazioni dove sono richieste alte prestazioni.

E' provvista di un'interfaccia USB2.0 ad alta velocità e di un ampio sensore CMOS che garantisce un ampio campo di vista.

Tali caratteristiche, assieme al software incluso BEL Eurisko per l'analisi delle immagini, permettono di catturare di eccellente qualità, con ottima resa di colori e risoluzione.

DV-5000	
Sensor - Sensore	1/2.5" color CMOS 5.0MPixel
Pixel size - Dim.Pixel	2.2 µm × 2.2 µm
Imaging area - Area immagine	5.7mm x 4.28mm
Sensitivity - Sensibilità	0.5 V/lux-sec
Resolution - Risoluzione	2592 × 1944
Speed - Velocità	3 f/s @ 2592×1944 12 f/s @ half resolution
White balance Bilanciamento bianco	Auto / Manual
Exposure control Controllo esposizione	Auto / Manual
Output - Uscita	USB2.0, 480Mb/s
Power supply - Alimentazione	USB2.0
Working Temperature - Temperatura d'esercizio	0°C ~ 50°C
Mount - Attacco	C mount

Package Includes:

Camera, USB cable, 1x C-mount adapter, Adapter 0.5x for eyepieces (biological and stereo), calibration micrometric slide, BEL Eurisko Image analysis software.

Incluso nella confezione:

Telecamera, Cavo USB, adattatore 1x di tipo C, Adattatore 0.5x per oculari (biologico e stereo), vetrino micrometrico di calibrazione, Software di analisi immagini BEL Eurisko.

Accessories

CCTV objective for environment shoots, cleaning kit

Accessori

Obiettivo CCTV per riprese d'ambiente, kit di pulizia

Ordering codes - Codici d'ordine	
DV cameras - Telecamere DV	
DV-5000	5.0MegaPixel, USB2.0. Windows XP/VISTA, software BEL Eurisko included <i>Telecamera 5Mpixel, USB2.0, Windows XP/VISTA, software BEL Eurisko compreso</i>
MB060025	CCTV objective for environment shoots <i>Obiettivo CCTV per riprese d'ambiente</i>
MB120010	Cleaning kit [100 tissue+solvent] <i>Kit pulizia [100 salviette+solvente]</i>

